

Needles & Syringes

Max-I-Probe-Dentsply MPL

Endodontic/Periodontal Irrigation Probes. Fully irrigate root canals (including the apical third). Closed, rounded tip and side dispersal port.

Endodontic Irrigation Kit

Contains: 5 Disposable Syringes 2 (3cc) & 3(5cc) and 20 Probes 5 (25g), 5 (28g), & 10 (30g)
MAXP005 55.35

Max-I-Probe Endodontic Irrigation Probe-Dentsply MPL

The cannula's closed-end, side window dispersal point allows safe cleaning of the gingival sulci.

- 100/Package
- Each 107.45
- MAXP211 Grn, 21G x 1" 107.45
- MAXP231 Lt. Bl, 23G x 1" 107.45
- MAXP241 Purple, 24G x 1" 107.45
- MAXP251 Org, 25G x 1" 107.45
- MAXP281 Red, 28G x 1" 107.45
- MAXP301 Dk Bl, 30G x 1" 107.45

Endo/Perio Starter Kit-Dentsply Rinn

Kit
Includes: 30 Probes
MAXP001 55.35

30 GA Endodontic Irrigation Kit
Includes: 30-30GA Probes
MAXP006 55.35

Endodontic Irrigation Needle-Kendall

Rigid polypropylene hard pack, tamper-evident heat stake, flow control notched tip, and epoxy needle attachment.

- 23 GA-1.25-31.7mm-Yellow
- 25/Box
- 8881-471232 1-8.95 4-8.55
- 27 GA-1.25-31.7mm-Yellow
- 25/Box
- 8881-471273 1-8.95 4-8.55

Sterile Syringes-Kendall

- 3cc-Luer Lock Syringes**
- 100/Box
- 8881-513934 1-12.85 10-11.65
- 6cc-Regular Luer Lock Syringes**
- 50/Box
- 8881-516937 1-12.85 10-11.65

Luer Lock Syringe with Endodontic Needle-Kendall

Sterile, Luer-Lock needle/syringe combo, disposable, with a length of 1 1/4", 31.7mm.

- 23 GA-Orange
- 100/Box
- 8881-513843 40.95
- 27 GA-Yellow
- 100/Box
- 8881-513850 40.95

Blunt Perio Irrigation Needles-Kendall

Sterile, aluminum Luer-Lock hub with autoclavable hard pack.

- 23 GA-1"-25.4mm-0.25"-Orange
- 100/Box
- 8881-202397 18.95

Pulp Canal Sealers & Cements

GuttaFlow-Coltene Whaledent

A new self-curing obturation material for root canals that combines two products in one: Gutta percha and sealer. Requires no heat.



Intro Kit
Contains: 20 GuttaFlow Capsules; 20 Canal Tips; 1 Dispenser; 100 Gutta Percha Point; Iso Assorted; #15-40; 60 Gutta Percha Points; .04 Greater Taper, Assorted #20-45
H08035 275.95

Refill Kit
Contains: 20 GuttaFlow Capsules; 20 Canal Tips
H08085 104.55
GuttaFlow Dispenser
H08008 157.45

AH-26R-Dentsply Maillefer

Used as a root canal sealer in conjunction with gutta percha, silver points or a solid core obturator.

Complete Package
667000 94.55

AH Plus Dentsply Maillefer

Complete Package
667003 103.35

Seal & Protect-Dentsply Maillefer

Sealant for exposed dentin offers the dentist a long-term solution. May be applied every 6 months. Provides a polymerized layer mechanically reinforcing the dentin to protect against erosion.

658010 73.85

Lentulo® Spiral Fillers-Dentsply Maillefer

Has consistently spaced spirals which provide outstanding flexibility in distributing root canal sealer and cement evenly throughout the root canal system.

- 25mm-4/Package
- Each 1-11.45 4-10.95
- 671570 Red #1
- 671571 Blue #2
- 671572 Green #3
- 671573 Black #4
- 671574 Assorted

- 21mm-4/Package
- Each 1-11.45 4-10.95
- 671560 Red #1
- 671561 Blue #2
- 671562 Green #3
- 671563 Black #4
- 671564 Assorted

- 17mm-4/Package
- Each 1-11.45 4-10.95
- 671550 Red #1
- 671551 Blue #2
- 671552 Green #3
- 671553 Black #4
- 671554 Assorted

Endodontic Sealers

WE SEARCH FOR LOW PRICES SO YOU DON'T HAVE TO!



Endodontic Sealers

T.E.R.M Temporary Endo Restorative-Dentsply Maillefer

Light-cured, resin based temporary restorative.

Complete Kit

Contains: 30 Compules® Tips-0.25g EA; 1 Compulex® Syringe.

670300 44.75

Refill Kit

Contains: 30 Compules® Tips-0.25g EA

670302 31.85

EZ-Fill-Essential Dental

The E-Z-Fill bi-directional spiral and epoxy root canal cement combined with a single point technique has been shown to create a seal equivalent to lateral condensation and thermoplastic gutta percha.

INTRO KITS CONTAINS:

4-Bi-Directional Spirals; 7.5gm Root Canal Epoxy Cement Gel; 8.0gm Powder; Measuring Scoop

Stainless Steel Intro Kit

1600-00 91.35

Nickel Titanium Intro Kit

1605-00 99.95

Refill Kit Stainless Steel

Contains: 1-21mm; 3-25mm length

1600-01 21.25

Refill Kit Nickel Titanium

Contains: 1-21mm; 3-25mm length

1605-01 29.15

EZ-Fill Epoxy Root Canal Cement

1608-00 74.65

EZ-Fill Gel Tube (ONLY) 7.5g

1610-00 32.65

CRCS® Calciobiotic Root Canal Sealer-Hygenic

CRCS® sealer, with calcium hydroxide, fulfills all the requirements of the "ideal sealer" and maintains its superiority long term.

36 Unit Complete Package

H03534 46.55

Sealapex (Non-Eugenol)-Sybron Endo

Radiopaque and easy to use. 45 minutes setting time after placement. Minimum of 2 hours working time.

Standard Package

Contains: Base (12gm); Catalyst (12gm)

18432 45.75

Pulp Canal Sealer-Sybron Endo

Working time of 20-40 minutes, original zinc oxide formula, radiopaque, and low solubility in oral fluids.

Standard Package

Contains: 10.5gm Powder; 4.0ml Liquid

24875 133.45

Pulp Canal Sealer EWT-Sybron Endo

Extended working time of 6-8 hours, radiopaque, and low solubility in oral fluids.

Standard Package

Contains: 10.5gm Powder; 4.0ml Liquid

24746 136.65

Tubli-Seal®-Sybron Endo

A two paste system used with Kerr Silver or Gutta Percha points for filling and sealing root canals.

Standard Package

Contains: 1 (10gm) Tube Base; 1 (3.5gm) Tube Accelerator; 1 Mixing Pad

00340 46.55

Tubli-Seal® EWT-Sybron Endo

Standard Package

25903 46.25

Multi-Cal-Pulpdent

All purpose calcium hydroxide preparation.

Kit

Contains: 4x (1.2ml) Syringes & 8 Applicator Tips

MULTI 31.35



Pulp Capping Paste-Pulpdent

For vital pulp and root canal therapy.



Pulpdent Paste Kit

Contains: 3ml Syringe with 24 Needles (18 GA Blunt)

PSYK 29.95

TempCanal-Pulpdent

Temporary calcium hydroxide root canal treatment paste.



Kit

Contains: 3ml Syringe & 12 Needles (22 & 25GA)

PTCK 31.95

CAVIT®-3M ESPE

A temporary filling material which hermetically seals in medicament.

4 Tubes

44016 17.45

20 Tubes

44024 77.85

28gm Jar

44030 17.45

Cavit G-28gm Jar

44313 17.45

RC-Prep™-Premier

Provides rapid cleansing and enlarging of root canals.



Jar 18gm

9007131 28.85

RC Prep Box

2 Syringes 9gms EA

9007135 30.35

RC Prep Syringe Kit

5-3cc Syringes & 50 Tips

9007129 34.65

Tip Refills-50 Angled Tips

9007139 15.45

Endodontic Sealers

Dental Supply and Equipment Can Help You Design or Remodel Your Office.

Savings!!

RC-Prep™ MicroDose™-Premier

Cleans and lubricates for chemo-mechanical preparation of root canal.

Intro Kit

Contains: 100 Pipettes and 2 Autoclavable Syringes.

9007138 94.95

Refill Kit

Contains: 100 Pipettes

9007137 80.95

Syringe Package

Contains: 3 Autoclavable Syringes

4009099 28.15

U/P Root Canal Cement-Sultan Chemists

Radiopaque, non-staining, non-irritating; 20 minute working time. Premeasured capsules.

Small Bulk Set

1/2 oz. Powder-1/4oz Liquid

11615 79.45

Endodontic Solutions, Prepackaged Disposable

Endodilator Juniors-Miltex

Contains: 15% EDTA

Pipettes(Green)-35/Box

017-39211P 47.55

Endo-Klean Juniors-Miltex

Contains: 2.5% Sodium Hypochlorite

Pipettes (Red)-35/Box

033-39214P 35.95

Chlorhexidine Juniors-Miltex

Contains: 2% Chlorhexidine

Pipettes(Lt Blue)-35/Box

033-39212P 49.95

Slide Juniors-Miltex

Contains: File Lubricant with 10% urea peroxide

Pipettes(Clear)-35/Box

033-39210P 47.55

Slide SH Juniors-Miltex

Contains: Slide with 3% Sodium Hypochlorite

Pipettes(Blue)-35/Box

017-39211P 47.55

Lube Juniors-Miltex

Contains: File lubricant w/NaOH & EDTA

Pipettes(Orange)-35/Box

033-39215P 54.55

Endo Medicaments

Glyde File Prep-Dentsply Maillefer

Syringe Kit

Contains: 5 (3gm) Syringes and 100 Disposable Tips

673021 32.55

Tip Refill

Contains: 100 Disposable Syringe Tips

673022 16.95

Camphorated Parachlorophenol-EPR

1oz Bottle

10302 15.45

Eugenol (USP)-EPR

1oz Bottle

10502 4.55

4oz Bottle

10504 13.75

FormoCresol-EPR

1oz Bottle

10601 15.15

Superoxol-EPR

Bleaching Product. 1oz

11513 13.85

EDTA Solution-Vista 17% Solution.

Chelating agent assists in instrumentation of the root canal.

120ml Bottle

317001 10.75

CHX Endo Solution-Vista

A 2%Chlorhexidene Gluconate Solution used to final irrigation step for long lasting disinfection.

Kit

Contains: 4 - 12ml prefilled syringe; 20-27ga irrigating tips.

502400 9.75

CHX Endo Refill

503350 15.25

Endosolv E Root Canal Filling Remover-Septodont

Softens filling materials for easy removal and improving access to root canal.

.44 fl. oz. Bottle

01-S0260 20.35

Chloroform-Sultan Healthcare

8 oz. bottle

20390 51.65

Endomet Root Canal Sealer-Septodont

Acts as an antiseptic and protects against infection.

Kit

Contains: Endomet Liquid (0.33 fl. oz); Endomet Powder (0.49 oz)

01-S0250 65.85

Buckley's Formo Cresol-Sultan Chemists

1oz

10203 37.25

Calcium Hydroxide Powder-Sultan Chemists

2 oz. Powder

10301 13.45

Endo-Ice®-Hygenic

Pulp vitality refrigerant spray-6oz Can

H05032 18.45

SmearClear-Sybron Endo

A 17% EDTA solution.

Kit

Contains: Solution; 3 Empty Syringes and 30 Side-Vented Needles.



973-1001 35.25

Bottle Refill

973-1002 29.95

3 Empty Syringes and 30 Needles

973-1003 8.95

READY TO PLACE YOUR ORDER?
HAVE ANY QUESTIONS?

Endodontic Sealers