

**Operative Pins with Accessories**

**TMS Thread Mate System®-Coltene Whaledent**

Gold plated stainless steel. Manual placement with hand driver or mechanical placement with AutoKlutch® chucks.

**TMS Regular .027" Gold**

Contains: 20 Self-Threading Pins, 1 Universal Hand Driver and 1 Kodex® Drill K-96

**T-01 ..... 103.75**

**T-02 Refill Kit: 50 Self-Threading Pins**

**T-02 ..... 163.85**

**T-03 Bulk Kit: 100 Self-Threading Pins**

**T-03 ..... 290.95**

**TMS Minim® .021" Silver**

Contains: 20 Self-Threading Pins, 1 Universal Hand Driver, and 1 Kodex® Drill K-92

**T-04 ..... 103.15**

**T-05 Refill Kit: 50 Self-Threading Pins**

**T-05 ..... 163.85**

**T-06 Bulk Kit: 100 Self-Threading Pins**

**T-06 ..... 290.95**

**TMS Self-Shearing (Double Shear)**

**TMS Regular™ .027" Gold**

T-08 Refill Kit: 100 Single Pins

**T-08 ..... 237.85**

**TMS Minim® .021" Silver**

T-11 Refill Kit: 100 Single Pins

**T-11 ..... 235.85**

**TMS Link Series®-Coltene Whaledent**

Self-threading and self-shearing pins used for the retention of restorations and build-up techniques.

**TMS Link® Minuta-(.0135")-Pink**

Complete Kit Contains: 20 Single Shear Pins, 1 Universal Hand Driver and 1 Kodex® Drill K-89

**L-501 ..... 123.35**

**TMS Link® Minikin-(.017")-Red**

Complete Kit Contains: 20 Single Shear Pins, 1 Universal Hand Driver and 1 Kodex® Drill K-91

**L-511 ..... 123.35**

**Minikin-(.017")-Red**

Bulk Kit Contains: 50 Self-Threading Pins

**L512 ..... 247.55**

**TMS Link® Minim-(.021")-Silver**

Complete Kit Contains: 20 Single Shear Pins, 1 Universal Hand Driver, and 1 Kodex® Drill K-92

**L-521 ..... 124.75**

**Minim-(.021")-Silver**

Bulk Kit Contains: 50 Self-Threading Pins

**L-522 ..... 184.75**

**TMS Link® Minim-(.021")-Silver-Double**

Complete Kit Contains: 20 Self-Threading Pins, 40 Single Pins, 2 Kodex® Drills K-92 and 1 Universal Hand Driver

**L-531 ..... 140.95**

**L-532 Bulk Kit**

Bulk Kit Contains: 50 Self-Threading Pins and 100 Single Pins

**L-532 ..... 242.95**

**TMS Link® Regular-(.027")-Gold**

Complete Kit Contains: 20 Self-Threading Pins, 1 Universal hand Driver, and 1 Kodex® Drill K-96

**L-541 ..... 123.95**

**TMS Link® Minim-(.027")-Gold-Double**

Complete Kit Contains: 20 Self-Threading Pins, 40 Single Pins, Universal Hand Driver, and 2 Kodex® Drills K-96

**L-551 ..... 140.95**

**TMS link® Minikin-(.017")-Red**

Refill Kit: 50 Single Shear Pins

**L-812 ..... 244.15**

**TMS Link Plus Series®-Coltene Whaledent**

**TMS Link Plus® Minim-(.021")-Silver**

Complete Kit Contains: 2 Single Shear Pins, 1 Universal Hand Driver and 1 Kodex® Drill K-92

**L-721 ..... 106.75**

**L-722 (Refill 50).....186.50**

**L-731 Refill**

Refill Kit Contains: 50 Single Shear Pins

**L-731 ..... 142.75**

**TMS Link Plus® Minim 2-in-1 Silver**

Refill Kit: 50 Double Shear Pins

**L-732 ..... 248.95**

**TMS Link Pluz® Regular-(.027")-Gold**

Bulk Kit Contains: 20 Single Shear Pins, 1 Universal Hand Driver, and 1 Kodex® Drill Kit

**L-741 ..... 106.75**

**TMS Link Plus® Regular 2-in-1 Gold**

Bulk Kit Contains: 20 Double Shear Pins, 1 Universal Hand Driver and 1 Kodex® Drill K-96

**L-751 ..... 145.45**

**L-752 Refill Kit**

Refill Kit Contains: 50 Double Shear Pins

**L-752 ..... 248.95**

**TMS Link Plus® Minim 2-in-1 Silver**

Bulk Kit Contains: 20 Double Shear Pins, 1 Universal Hand Driver, and 1 Kodex® Drill K-92

**L-841 ..... 14.35**

**Max® Restorative Pins-Coltene Whaledent**

Provides increased retention of restorations and build-up techniques.

**Max™017 Blue**

Complete Kit Contains: 20 Max™ Restorative Pins, 2 Max™ Drills M-19, 1 Universal Hand Driver and 1 Steri-Cord®

**M-11 ..... 124.75**

**Max™017 Blue**

Bulk Kit Contains: 50 Max™ Restorative Pins, 5 Max™ Drills M-19, 1 Universal Hand Driver, and 1 Steri-Cord®

**M-13 ..... 185.95**

**Max™021 Purple**

Complete Kit Contains: 20 Max™ Restorative Pins, 5 Max™ Drills M-92, 1 Universal Hand Driver, and 1 Steri-Cord®

**M-21 ..... 122.75**

**Max™021 Purple**

Bulk Kit Contains: 50 Max™ Restorative Pins, 5 Max™ Drills M-92, 1 Universal Hand Driver, and 1 Steri-Cord®

**M-23 ..... 185.95**

**Filpin™ System-Filhol Dental**

Pure Titanium Self-Shearing retention pins for automatic and manual placement. Filpin Standard Package: 35-Pins .60mm/.021 Universal, 1 Filpin Universal Twist Drill.

**Standard Package Universal**

**PIN101 ..... 94.75**

**Standard Package Large**

**PIN102 ..... 94.75**

**Endo Posts**

**Flexi-Post®-Essential Dental Systems**

The split shank closes during insertion, absorbing insertional stress.

**Assorted Kit 0-1-2**

**110-00 ..... 183.85**

**Assorted Kit 1-2-3**

**110-01 ..... 183.85**

## Pins & Posts

### Flexi-Post®-Essential

#### Refills

<b>Size 00-White</b>	
130-00 .....	125.95
<b>Size 0-Yellow</b>	
130-0 .....	125.95
<b>Size 1-Red</b>	
130-01 .....	125.95
<b>Size 2-Blue</b>	
130-02 .....	125.95
<b>Size 3-Green</b>	
130-03 .....	125.95

### ParaPost Fiber White-Coltene Whaledent

A great alternative to metal posts when a metal-free restoration is desired.

#### Intro Kit

Contains: 10 Assorted posts, 4 Drills, 1 Hand Driver

PF-160 ..... 150.45

#### Refill Kits-Qty. 5 Pieces

Each .....	45.55
PF161-4.5 .....	Size 4.5
PF161-5 .....	Size 5
PF161-5.5 .....	Size 5.5
PF161-6 .....	Size 6

### IntegraPost™ System-Premier

A secure stress-free antirotational core.

#### Intro Kit

Contains: 18 Posts, 3 EA #3-#5, 5 EA #4-#4.5-, 2 EA #6, 5 Drills Sizes 1-6, IP Carrier, IP Probe, and Autoclavable Hot Box

3001002.....272.95

#### Mini Intro Kit

Contains: 6 Posts: 3 EA #3-#4-#4.5, 3 Drills EA #3-#4-#4.5, IP Carrier and Lucite Block

3001001 ..... 86.75

#### IntegraPost Refills-10/Package

Each .....	81.55
3001105 .....	Size 3
3001110 .....	Size 4
3001115 .....	Size 4.5
3001120 .....	Size 4.5L
3001125 .....	Size 5
3001130 .....	Size 6



### Integra Drills-3/Package

Each .....	40.35
3001205 .....	Size 3
3001210 .....	Size 4
3001215 .....	Size 4.5
3001225 .....	Size 5
3001230 .....	Size 6

### Assorted-5/Package

3001235 .....	58.45
<b>IP Carrier</b>	
3001301 .....	13.45
<b>IP Probe</b>	
3001305 .....	21.35

### Flexi-Flange®-Essential Dental Systems

#### Stainless Steel Intro Kit

Contains: 12 Serrated Posts w/Flange (4 Posts of 3 Sizes), 3 Primary Reamers, 3 Countersink Drills, 1 Extender, 1 Wrench, 1 Depth Gauge, 5 Gates Glidden Drills

410-00 ..... 186.85

#### Stainless Steel Refill Package

Contains: 10 Serrated Posts w/Flange, 1 Primary Reamer, 1 Countersink Drill

Each .....	128.95
430-0 .....	Yellow/Size 0
430-01 .....	Red/Size 1
430-02 .....	Blue/Size 2
430-03 .....	Green/Size 3

### ParaPost Fiber

#### Lux-Coltene Whaledent

A translucent fiber resin material improves esthetics for anterior teeth, minimizes shadowing. Ability to transmit polymerization light allows for the convenience of using self, dual or light cured materials.



**Dental Equipment and Supply Can Install New Chairs, Cabinets, Lights and Much More!**

### Intro Kit

Contains: 2 Posts (Size 3), 1 (#3) Drill, 3 Posts (Size 4), 1 (#4) Drill, 3 Posts (Size 4.5), 1 (#4.5) Drill, 3 Posts (Size 5), 1 (#5) Drill, 2 Posts (Size 5.5), 1 (#5.5) Drill 2 Posts (Size 6), 1 (#6) Drill

PF-170 ..... 193.75

#### Refills

Each .....	45.85
PF171-3 .....	Size 3
PF171-4 .....	Size 4
PF171-4.5 .....	Size 4.5
PF171-5 .....	Size 5
PF171-5.5 .....	Size 5.5
PF171-6 .....	Size 6

### Para-Post XT-Coltene Whaledent

Permits shortening without compromising mechanical retention.



#### Intro Kit

Contains: 3 Posts (Size 3), 5 Posts (Size 4), 5 Posts (Size 4.5), 4 Posts (Size 5.5), 3 Posts (Size 6), 6 Drills 1 of EA Side, 1 Drill Driver, 1 Post Driver (For 3, 4, 4.5 Sizes), 1 Post Driver (For 5, 5.5, 6 Sizes)

P-680T ..... 302.95

#### Refills

<b>.036"-Brown-10 Pieces</b>	
P-683-0 .....	93.45
<b>.040"-Yellow-10 Pieces</b>	
P-684-0 .....	93.95
<b>.045"-Blue-10 Pieces</b>	
P-684-5 .....	93.95
<b>.050"-Red-10 Pieces</b>	
P-685-0 .....	93.95
<b>.055"-Purple-10 Pieces</b>	
P-685-5 .....	93.95
<b>.060"-Black-10 Pieces</b>	
P-686-0 .....	93.95

### ParaPost® XP™-Coltene Whaledent

One-Office-Visit Intro Kit. Improved Para Post design for amalgam cores and multi-rooted teeth.



#### Intro Kits

<b>Stainless Steel</b>	
P-780 .....	248.55
<b>Titanium Alloy</b>	
P-780T .....	229.55

**Stainless Steel Refill Kits**

.036"-Brown-10 Pieces P-744-3 .....	73.95
.040"-Yellow-10 Pieces P-744-4 .....	73.95
.045"-Blue-10 Pieces P-744-4.5 .....	73.95
.050"-Red-10 Pieces P-744-5 .....	73.95
.055"-Purple-10 Pieces P-744-5.5 .....	73.95
.060"-Black-10 Pieces P-744-6 .....	73.25
.070"-Green-10 Pieces P-744-7 .....	73.25

**Titanium Alloy Refill Kits**

.036"-Brown-10 Pieces P-784-3 .....	74.95
.040"-Yellow-10 Pieces P-784-4 .....	74.95
.045"-Blue-10 Pieces P-784-4.5 .....	74.55
.050"-Red-10 Pieces P-784-5 .....	74.95
.055"-Purple-10 Pieces P-784-5.5 .....	74.95
.060"-Black-10 Pieces P-784-6 .....	74.95
.070"-Green-10 Pieces P-784-7 .....	74.95

**ParaPost® XPT™ Casting Technique Kit-Coltene Whaledent**

<b>Intro Kit</b> 25 EA ParaPost XP Impression: Temporary, Burnout Posts, 7 ParaPost X Drills-1EA P-781 .....	251.55
--	--------



**ParaPost®XH™  
Coltene  
Whaledent**

One-Office-Visit  
Intro Kit. Titanium Alloy.



<b>Intro Kit</b> P-880 .....	276.45
---------------------------------	--------

**Refill Kits**

.036"-Brown-10 Pieces P-88-3 .....	79.85
.040"-Yellow-10 Pieces P-88-4 .....	79.85
.045"-Blue-10 Pieces P-88-4.5 .....	79.85
.050"-Red-10 Pieces P-88-5 .....	79.85
.055"-Purple-10 Pieces P-88-5.5 .....	79.85
.060"-Black-10 Pieces P-88-6 .....	79.85
.070"-Green-10 Pieces P-88-7 .....	79.85

**ParaPost® System-  
Coltene Whaledent**

Universal endodontic post system for direct post and core build-ups and all casting techniques. Optimum of functional forces.

**Drills**

.036"-Brown-3 Pieces P-42-3 .....	52.55
.040"-Yellow-3 Pieces P-42-4 .....	52.55
.050"-Red-3 Pieces P-42-5 .....	52.55
.060"-Black-3 Pieces P-42-6 .....	52.55
.070"-Green-3 Pieces P-42-7 .....	52.55



**Plastic Impression Posts**

.036"-Brown-3 Pieces P-43-3 .....	35.25
.040"-Yellow-3 Pieces P-43-4 .....	35.25
.050"-Red-3 Pieces P-43-5 .....	37.25
.060"-Black-3 Pieces P-43-6 .....	37.25

**Para-Post® Accessories-  
Coltene Whaledent**

<b>Universal Hand Driver (1)</b> L-590 .....	5.15
<b>Metal Wrench A</b> For Sizes 3, 4, 4.5 Posts L-594 .....	30.95
<b>Metal Wrench B</b> For Sizes 5, 5.5, 6 Posts L-595 .....	30.95

**Pin Drills**

**Kodex Drills - Coltene Whaledent**

Manufactured from high-speed tool steel and designed for use with all Whaledent self threading pin systems.or Sizes 5, 5.5, 6 Posts  
6/pk .....

K90.....	.017" x .059: SS Red
K91.....	.017" x .059: LS Red
K92.....	.021" x .079: Silver
K93.....	.021" x .157: Silver
K96.....	.027" x .079: Gold
K97.....	.027" x .197: Gold
K98.....	.028" x .197: Green



**DON'T SEE WHAT  
YOU ARE  
LOOKING FOR?**

**LET US HELP....CALL  
NOW  
319-362-2265  
OR  
1-800-365-6265**

Pins & Posts