

**Retraction Materials**

**Stay-Put Retraction Cord-  
Coltene Whaledent**

Unique combination of softly braided retraction cord and an ultrafine copper filament that holds its shape. Non-impregnated.



Each .....	13.55
RE521000 .....	0 X-Fine
RE521001 .....	1 Fine
RE521002 .....	2 Medium
RE521003 .....	3 Thick

**Z-Twist™100% Cotton Weave-  
Gingi-Pak**

Aluminum Chloride-Teal-72"

Each .....	1-8.95	4-7.85
12170 .....	00 Very Thin	
12171 .....	1 Thin	
12172 .....	2 Medium	
12173 .....	3 Thick	

Non-Impregnated Purple

Each .....	1-7.95	4-6.95
11170 .....	00 Very Thin	
11171 .....	1 Thin	
11172 .....	2 Medium	
11173 .....	3 Thick	

Epinephrine-Green

Each .....	1-8.95	4-7.95
10170 .....	00 Very Thin	
10171 .....	1 Thin	
10172 .....	2 Medium	
10173 .....	3 Thick	

**Original Cord-Gingi-Pak**

**Epinephrine: 2 Ply.**

10120 .....	1-8.95	4-7.95
-------------	--------	--------

**Gingi-Plain: 2 Ply**

11120 .....	1-8.85	4-7.95
-------------	--------	--------

**Soft Twist Cord-Gingi-Pak**

Epinephrine: 72"

Each .....	1-8.85	4-8.45
10105 .....	1 Thin	
10110 .....	2 Medium	
10115 .....	3 Thick	

Non-Impregnated

Each .....	1-7.85	4-7.05
11105 .....	1 Thin	
11110 .....	2 Medium	
11115 .....	3 Thick	

**Gingi-Pak Pellets with Epinephrine-  
Gingi-Pak**

**Package**

10135 .....	17.65
-------------	-------

**Cotton Coil w/Epinephrine**

10140 .....	17.65
-------------	-------

**Standard Cord Packing Instruments-  
Pascal**

A unique shoulder curve at both ends.

**Plain: R-101**

26-100 .....	20.55
--------------	-------

**Serrated: R-11**

26-110 .....	20.55
--------------	-------

**#113 Gingival Cord Packer-  
Hu-Friedy**

GCP113 .....	25.45
--------------	-------

**CSII Gingival Cord Packer-  
Hu-Friedy**

GCPCSI1 .....	25.45
---------------	-------

**EXPA-SYL™ Gingival Retraction  
System-Kerr**

A fast, safe, and painless alternative to packing retraction cords. Expa-syl is a paste for gingival retraction that not only opens the sulcus but also leaves the field dry.



**Starter Kit**

Contains: 6 Capsules Retraction Past (1g EA), 12 Applicator Tips, 1 Applicator Gun

32379 .....	337.45
-------------	--------

**Capsules-20/Package**

31172 .....	289.95
-------------	--------

**Applicator Tips-40/Package**

31173 .....	81.65
-------------	-------

**Applicator Tips-100/Package**

32264 .....	147.25
-------------	--------

**Applicator Gun**

31174 .....	319.65
-------------	--------

**Gel-Cord™-Pascal**

Aids in gingival retraction and controls bleeding. Unit dose cartridges of 25% Aluminum sulfate gel with a raspberry flavor.

**Intro Kit**

Contains: 25 (0.32gm) Cartridges, 25 Applicator Brushes, Bottle Plain Cord

15-500 .....	79.95
--------------	-------

**Refill Kit**

Contains: 25 (0.32) Cartridges

15-505 .....	56.95
--------------	-------

**Bulk Kit**

Contains: 3 (10ml) Syringes, 25 Empty Tubes, 25 Applicator Brushes

15-525 .....	128.25
--------------	--------

**Pro Pack 12 Syringe**

Contains: 12 (1g) Syringes and 14 Tips

15-540 .....	37.95
--------------	-------

**Pro Pack 40 Syringe**

Contains: 40 (1g) Syringes and 50 Tips

15-540 .....	102.85
--------------	--------

**Racord-Pascal**

Twisted, medicated, retraction cord.

Each .....	9.75
07-700 .....	8
07-710 .....	9
07-720 .....	10

**Hemodent™-Premier**

Buffered aluminum chloride, epinephrine-free liquid.



**10cc Bottle**

9007071 .....	1-14.55	4-13.99
---------------	---------	---------

**20cc Bottle**

9007072 .....	24.95
---------------	-------

**Hemodent Cord-Premier**

Epinephrine free hemodent cords are already impregnated with hemodent liquid.



Easily packed subgingivally without fraying or disintegrating.

Each .....	13.55
9007077 .....	Regular Twisted #3
9007078 .....	Regular Twisted #9
9007065 .....	Braided Thin
9007067 .....	Braided Med/Thin

**Hemogin-L-Van R**

25% aluminum chloride in a lubricating glycerin base. Hemogin-L, 10cc bottle w/eye dropper top.



13057 .....	15.25
-------------	-------

**Hemogin-L: 40cc Bottle**

13052 .....	41.95
-------------	-------

**CALL NOW !!!!!**

Retraction Materials

## Retraction Materials

### GingiGel-DUX

The braid is packaged in and pulled through a gel containing 20% buffered aluminum chloride. 10ft. of braid per bottle.



Each ..... 13.25  
 13165 ..... 1 Small Size  
 13166 ..... 2 Medium Size  
 13167 ..... 3 Large Size

### GingiBraid-DUX

Braided retraction materials that will not split or collapse in the sulcus. Epinephrine/Alum 87  
 Gingibraid/A solution of approximately 8% Racemic Epinephrin (a vasoconstrictor) and 7% Aluminum Potassium Sulfate (an astringent), impregnated and dried. 72 inches.



Each ..... 12.95  
 13170 ..... 0 Fine Size-Blue  
 13171 ..... 1 Sm. Blue w/Purple Strand  
 13172 ..... 2 Med. Blue w/Green Strand  
 13173 ..... 3 Lg. Blue W/Pink Strand

### Aluminum Potassium Sulfate GingiBraid

A 10% solution, impregnated and dried.

Each ..... 12.95  
 13210 ..... 0 Fine Size-Yellow  
 13211 ..... 1 Sm. Yellow w/Purple Strand  
 13212 ..... 2 Med. Yellow w/Green Strand  
 13213 ..... 3 Lg. Yellow w/Pink Strand

### Non-Impregnated GingiBraid

Each ..... 11.75  
 13217 ..... 0 Fine Size-White  
 13214 ..... 1 Sm. White w/Purple Strand  
 13215 ..... 2 Med. White w/Green Strand  
 13216 ..... 3 Lg. White w/Pink Strand

### Stasis®-Gingi-Pak

A 21% ferric sulfate hemostyptic formulation offering unsurpassed hemorrhage control, when measured against any non-epinephrine solution.

40ml  
 13115 ..... 32.95

### Styptin-DUX

20% buffered aluminum chloride hemostatic solution. A tasteless fast-acting hemostatic agent. The color (light blue) allows the operator to determine placement and from where it must be rinsed.



15ml Cranberry Styptin-Squeeze Bottle  
 13157 ..... 19.95

### GingiKnit-DUX

Constructed of thousands of tiny, absorbent loops knitted in long, interlocking chains. GingiKnite places easily without fraying or memory. 10 Feet of knit per bottle.



### Aluminum Potassium Sulfate GingiKnit

A 10% solution, impregnated and dried.

Each ..... 10.75

13509 ..... 00 Fine, Black/Orange/White  
 13510 ..... 0 Small, Orange/White  
 13511 ..... 1 Med. Black/Pink/White  
 13512 ..... 2 Lg. Pink/White

### Non-Impregnated GingiKnit

Each ..... 10.75  
 13499 ..... 00 Fine, Black/Orange/White  
 13500 ..... 0 Small, Orange/White  
 13501 ..... 1 Med., Black/Pink/White  
 13502 ..... 2 Lg., Pink/White

### Comprecap Compression Caps-Coltene Whaledent

Compression caps used to pen the sulcus, hold retraction cord in place and provide dry clean area for impression.



Assorted-Qty. 150  
 Contains: #1 (7mm), #3 (10mm), #5 (12.5mm)  
 RE530000 ..... 25.85

#1 Small-Qty. 120  
 RE530008 ..... 19.45

#3 Medium-Qty. 120  
 RE530010 ..... 23.65

#5 Large-Qty. 60  
 RE530014 ..... 16.65

#1 Small Bulk-Qty. 600  
 RE531008 ..... 81.85


#3 Medium Bulk-Qty. 600  
 RE531010 ..... 99.75

#5 Large Bulk-Qty. 400  
 RE531014 ..... 77.65

Retraction Materials

**MedConduit™**  
 An NDC Company

**WE NOW HAVE ONLINE ORDERING  
 COMPARE PRICES!! COMPARE PRODUCTS!!  
 SAVE TIME AND MONEY!!**



**GIVE US A CALL TODAY FOR MORE INFORMATION**