

Blades & Scalpels

Standard Handles-Miltex

Knife Handle #3

Fits surgical blades 10 thru 15. Extra Fine, Stainless.

4-7 13.25

Knife Handle #4

Fits surgical blades 20 thru 35. Extra Fine, Stainless.

4-8 13.25

Knife Handle #6

Fits surgical blades 20-25. Extra Fine, Red Plastic.

4-20 21.25

Scaples-Miltex

Disposable stainless steel. Sterile. 10/Box

Size 10-Standard

4-410 10.65

Size 11-Standard

4-411 10.65

Size 12-Standard

4-412 10.65

Size 15-Standard

4-415 10.65

Sterile Surgical Blades-Miltex

Disposable. Sterilized by Gamma Radition (2.5 M-Rads). Individually foil packaged, hermetically sealed blades.

Stainless Steel or Carbon Steel-100/Box

EachBox: 37.75

4-310 10

4-311 11

4-312 12

4-315 15

4-325 25

Special Blades

Stainless Steel or Carbon Steel-100/Box

Each 75.25

4-312B 12B

4-315C 15C

Protected Surgical Blades-Crosstex

50/box.....79.90

373910.....#10

373911.....#11

373915.....#15

Protected Surgical Blades-Crosstex

374030.....34.90

Quala Surgical Blade Handle #3



Stainless Steel

Q6001-3 3.75

Quala Surgical Blades

Stainless Steel



Sterile-100/Box

Each 24.95

Q3001-10 10

Q3001-11 11

Q3001-12 12

Q3001-15 15

Each 27.95

Q3001-12B 12B

Q3001-15C 15C

Quala Surgical Blades

Carbon Steel



Sterile-10/Box

Each 3.45

Q2001-10 10

Q2001-11 11

Q2001-12 12

Q2001-15 15

Each 5.25

Q2001-12B 12B

Q2001-15C 15C

Quala Disposable

Scalpel



10/Box

Each 7.15

Q6008-10 10

Q6008-11 11

Q6008-12 12

Q6008-15 15

Each 7.75

Q6008-12B 12B

Cold Packs

Instant Cold Pack-Kimberly Clark

4" x 10"-24/Case

59681 16.95

6.25" x 8.5"-24/Case

59688 16.95

Medicaments

ANTISEPTIC

Buckley's Formo-Cresol-Sultan Chemists

Contains 19% formaldehyde, 35% cresol, and 17.5% glycerin. Used to devitalize and sterilize infected pulp chambers and canals. 1oz.

10203 39.35

Camphorated Parchlorophenol (USP)-EPR-Sultan Chemists

Contains 35% parachlorophenol, 65% camphor. Used for disinfecting root canals and treating periapical infections. Virtually non-irritating. 1oz.



10302 15.45

Eugenol (USP)-EPR

Comparable to all leading brands.

1oz Bottle

10502 4.55

4oz Bottle

10504 13.75

Formo-Cresol-EPR

1oz Bottle

10601 15.15

BLEACHING

Superoxol-EPR

1oz

11513 13.85

DRY SOCKET

Surgical Products/Medicaments

**Sterilizer not quite right? Handpiece acting funny ?
Air compressor making strange noises ?
BEFORE YOU BUY CALL US!**

Surgical Products/Medicaments

Dry Socket Paste-Sultan Chemists

Used to prevent and treat boen and socket pain. Gives instant relief, washes out gradually as socket heals. 1oz.

10404.....
70.85



Alvogy-Septodont

Alveolar Dressing-11gm Jar
S0010 58.25

Dry Socket Remedy-Sultan Chemists

Liquid treatment introduced into socket on iodoform gauze or absorbable hemostatic material. 1oz bottle.

10405..... 71.55

HEMOSTAT

ActCel Hemostatic Gauze-ActSys Medical Inc.

The ActCel 2x2 Gauze is extremely versatile in its use and although the guaze can be used in any bleeding situation, it has been widely used in the dental industry to control bleeding under particular dental procedures. Dental Indications: Extraction sockets, Dry Socket Prevention, Periodontal Surgery, Endodontic Surgery, Other Oral surgery procedures.

HG020022 109.95

Surgicel Absorbable Hemostat-J&J

Assists in the control of capillary or venous bleeding or small arterial hemorrhage.

12 Envelopes-1/2 sq. in.
1955 274.35

Hemodettes Unit-Dose Cups-Dux

20% buffered aluminum chloride lubricating hemostatic gel packaged in unique plastic cups with attached finger ring. Unit-dose packaging eliminates cross-contamination. Light blue color indicates area of placement. The gel viscosity enhances placement and eliminates migration from the prep.



Cranberry Flavor-25/Jar
13150 19.95

Gel Foam Sterile Sponges-Pharmacia

Absorbable Gelatin Sponge-3/4" x 2 1/4" Strips-12/Box
PFZ1545 156.85

Gel Foam #4 Dental Sponge-Pharmacia

Absorbable Gelatin Sponge-1/2 sq. in.-Qty. 6
PFZ1575 74.15

Surgical Dressing & Periodontal Pack

Coe-Pak™ Automix NDS-GC America

Has no burning sensation, no unpleasant taste, no disagreeable odor, and offers proven protection to surgical sites.

Kit

Contains: 2 Double Barrel Cartridges (50ml EA), 12 Mixing Tips, and Instructions

135003 72.45

Bone Grafting

PerioGlas-Sunstar Butler

PerioGlas is the premier brand of bone graft material that aids in the re-generation of the bone and periodontal tissues, and has superior handling characteristics.



PerioGlas 6 x 0.5cc
18705P 280.55

PerioGlas 2 x 0.5cc
18721P 106.15

PerioGlas 1 x 1.5cc
18739P 129.35

BioPlant-Kerr

.50gm w/Biotip. Synthetic Bone. Patented slowly resorbable calcium hydroxide-coated copolymer. Promotes new bone growth. Prevents dry sockets.



.50gm Syringe Straight-Qty. 6
216104 670.95

.50gm Syringe Curved-Qty. 6
216114 670.95
.25 Syringe Curved-Qty.2

216132 131.95.
.25 Syringe Curved-Qty.6
216112 379.95

.125gm Syringe Curved-Qty.6
216110.....225.15

Bifoil Surgical Protective Strips-Kerr

6/Package
111012

58.85

Sutures

Quala Dental Sutures

The Quala labeled suture range covers the six most popular dental sutures. Box of 12.



18 Inch

3/0 SilkRev Cut: 1/2
Type: C26
P00062.....1-22.45 6-20.95

3/0 Silk Rev Cut: 1/2
Type: C31
P00064 1-21.45 6-20.75

3/0 Plain Gut Rev Cut: 1/2
Type: C31
P00065 1-24.55 6-23.75

4/0 Silk Rev Cut: 3/8
Type: C6
P00066 1-21.45 6-20.75

4/0 Plain Gut Rev Cut: 3/8
Type: C6
P00071 1-24.75 6-23.75

3/0 Silk Rev Cut: 3/8
Type: C6
P00072 1-21.45 6-20.65

27 Inch

4/0 Chromic Gut Rev Cut: 3/8
Type: C6
P00635 1-25.95 6-24.95

3/0 Chromic Gut Rev Cut: 3/8
Type: C6
P00636 1-25.95 6-25.95

3/0 Plain Gut Rev Cut: 3/8
Type: C7
P00810 1-25.95 6-25.95

3/0 Plain Gut Rev Cut: 3/8
Type: C6
P00066 1-25.95 6-25.95

Ethicon Sutures-Ethicon-J&J

Non Absorbable Silk Black Braided

4/0
Needle: P-3
Size: 18"
641G.....1-125.95.....6-122.55

4/0
Needle: C-3
Size: 18"
736G.....1-81.95.....6- 78.95

3/0
Needle: X-1
Size: 18"
632G.....1-68.85.....6-66.35

4/0
Needle: FS-2
Size: 18"
683G.....1-69.95.....6-66.65

5/0
Needle: FS-2
Size: 18"
682G.....1-72.25.....6-68.55

Absorbable Gut Medium Chromic

3/0
Needle: PS-2
Size: 18"
1638H.....3 Dozen.....484.95

4/0
Needle: PS-2
Size: 18"
1637G.....1 Dozen.....181.95

3/0
Needle: FS-2
Size: 27"
636H.....3 Dozen307.65

4/0
Needle: FS-2
Size: 27"
635H 3 Dozen 304.65

4/0
Needle: P-3
Size: 18"
1654G 164.65

5/0
Needle: P-3
Size: 18"
687G.....1-196.45.....6-191.95

Gut Plain

5/0
Needle: P-3
Size: 18"
686G.....1-195.45.....6-191.95

3/0
Needle: X-1
Size: 18"
612G.....1-89.35.....6-87.55

4/0
Needle: FS-2
Size: 27"
H822H 3 Dozen 268.65

4/0
Needle: FS-2
Size: 27"
H821H 3 Dozen 305.65

Coated Vicryl Violet Braided

3/0
Needle: FS-2
Size: 18"
J393H 236.95



4/0
Needle: FS-2
Size: 18"
J392H 236.95

4/0
Needle: P-3
Size: 18"
J464G 145.95

Perma Sharp™-Hu-Friedy

Stays sharp throughout procedures, holds shape, and have better point geometer.

4/0 Gut, Chromic
Needle: 3/8 Circle Premium Reverse Cut
Type: C3
PSN1654C 105.65

3/0 Silk, Black Braided
Needle: 3/8 Circle Reverse Cut
Type: C6
PSN7772S 44.25

4/0 Silk, Black Braided
Needle: 3/8 Circle Reverse Cut
Type: C6
PSN683S 44.55

4/0 Gut, Plain
Needle: 3/8 Circle Reverse Cut
Type: C6
PSN7771A 49.45

4/0 Gut, Plain
Needle: 3/8 Circle Reverse Cut
Type: C6
PSNH821A 49.45

3/0 Silk, Black Braided
Needle: 3/8 Circle Reverse Cut
Type: C7
PSN684S 44.55

3/0 Gut, Plain
Needle: 1/2 Circle Reverse Cut
Type: C9
PSN612A 50.15

3/0 Silk, Black Braided
Needle: 1/2 Circle Reverse Cut
Type: C9
PSN632S 44.55

3/0 Silk, Black Braided
Needle: 1/2 Circle Reverse Cut
Type: C43
PSN7762S 40.45

Surgical Products/Medicaments

SUTURE COMPARISON CHART

NEEDLE	QUALA	ETHICON	HU-FRIEDY
C26	P00062	2762	
C31	P00064	632G	
C31	P00065	612G	PSN612A
C6	P00066	683G	PSN683S
C6	P00071	7771	PSN7771A
C6	P00072	7772	PSN7772S
C6	P00635	635H	PSN635C
C6	P00636	636H	PSN636C
C7	P00810	H810H	PSNH810A
C6	P00822	H822H	PSNH822A

