

Waxes

Periphery Wax

Surgident® Periphery Wax-Heraeus Kulzer

For border molding. Specially formulated for hydrocolloid impressions. Tacky but not greasy.

Blue-60 Sticks/Box
50092189 19.95

Baseplate

Waxes-Dentsply Trubyte

Set-Up Wax

Dentsply® Set-Up Wax is specially formulated for arranging and articulating teeth. It carves easily and boils-out clearly.

Sheets-5 lb Package
7707901 59.35

Neowax Flexible Baseplate Wax

A popular all-season, all-purpose high grade synthetic baseplate wax that is tough, pliable and boils-out clearly.

0.050" Thick 5 lbs.
774140101 35.95

Truwax® Baseplate Wax

Truwax® is a high quality baseplate wax. It carves cleanly and sharply at room temperature and boils out clearly.

1 lb Package
Each 13.95
7709002 Regular
7708902 Extra Hard

5 lb Package
Each 49.55
7709003 Regular
7708903 Extra Hard

Medium Soft No. 3 Pink Wax-Hygenic

This wax is easily softened and cooled, meaning minimal shrinkage during set-ups. It takes a high polish allowing an accurate preview of the denture or partial denture case.

1 lb.
H00806 12.75
5 lb.
H00812 46.95

Modern Pink No. 3-Heraeus Kulzer



1 lb.
50093312 17.95
5 lb.
50093352 77.65

Shur Wax®-Heraeus Kulzer

For set-ups, occlusal rims, and final wax-ups. Extra tough, yet pliable all-purpose wax. Shaded to resemble gum tissue



Regular-1 lb.
50093152 13.45

X-Hard-1 lb.
50093252 13.45

Regular-5 lb.
50093112 63.45

X-Hard-5 lb.
50093252 63.45

Bite Waxes, Wafers, & Sheets

Denture Wax-Aluwax

Excellent for check bites. A box contains a minimum of 28 sheets.

DENTURE 21.85

Laminated Blue Bite Wafers-Hygenic

Aluminum foil laminated between two layers of soft bite wax lessens distortion and prevents patient from biting through and destroying the impression.

24/Box
H00825 14.75

Wax Bite Blocks-Heraeus Kulzer

Identical to high-quality Shur Wax® formulation. Pre-formed bite blocks save time. 11mm thickness.



Soft-Trial Package-12 Blocks
50095592 17.75

Soft-Lab Package-120 Blocks
50095992 61.95

Coprwx® Bite Wafers-Heraeus Kulzer

Thin foil layers prevent teeth from biting through. Copper particles provide uniform flow of heat and soften quickly.



84/Box
50092178 66.25

Sticky Waxes

Sticky Wax-Kerr

A tough wax for crown and bridge work. Also casting, soldering and a variety of other uses. A clean wax that leaves a residue-free burn-out.

Sticks-12/Small Box
00623 10.55

Sticks-120/Bulk Box
00625 35.35


Sticky Wax-Dentsply Trubyte

1 lb. Package of Sticks
77310 88.95

Boxing & Utility Waxes

Quala Boxing Wax

For boxing impressions. Soft and pliable. 1lb box. 1 1/2" x 11" sheets. Red.



Regular
50094698 13.75

X-Thin
50094699 13.75

Quala Utility Wax Ropes and Strips

For border molding impression trays and beading impressions to preserve the peripheral roll. Tacky but not greasy. Soft and pliable.

Each 1-10.95 4-9.95
50094955 Red Sq. Rope, 55/Box
50094954 White Sq. Rope, 55/Box
50094953 Lg. Red Strips, 55/Box
50094952 Lg. White Strips, 55/Box
50094951 Sm Red Strips, 114/Box
50094950 Sm White Strips, 114/Box

Waxes

MedConduit™
An NDC Company

**WE PASS
THE SAVINGS
ON TO YOU**

MedConduit™
An NDC Company